

Analysis suite

ProUmid GmbH & Co. KG
 August-Nagel-Straße 23 · 89079 Ulm · Germany
 Tel +49 (0)7305 2030 440 · Fax +49 (0)7305 2030 439
 info@proumid.com · www.proumid.com

proUmid GmbH & Co. KG


Optional software package for moisture sorption data import and advanced data analysis

Features


- ✓ Tabular and graphical display of data
- ✓ Configurable diagrams: sorption kinetics and isotherm
- ✓ Extended display options for multi-sample data
- ✓ Curve fitting: BET, GAB
- ✓ Calculation of the water vapor transmission rate of foils
- ✓ Tool for simplified evaluation of calibration measurements
- ✓ User configurable template-based report generator

License

- ✓ Data import from a single SPS/Vsorp system.
- ✓ The license is linked to the serial number of the device.
- ✓ Unlimited number of installations and users.



Clear presentation of measurement data



Easy curve fitting to common sorption models